

CSM-BC201 QUICK SETUP GUIDE

NOTE: Consult the "Directions For Use" included with monitor for complete instructions, troubleshooting, and safety information

Setup



- 1** Attach the Velcro strap.



- 2** Open battery compartment:
1. Install batteries
 2. Set delay switch to desired time (0, 2, or 4 seconds)
 3. If desired set "Pause" switch to "ON" (pressing "RESET" button allows patient transfer for 30 seconds without activating alarm).
 4. Close the cover and secure the screw



- 3** Write in the appropriate date(s) on sensor pad.



- 4** **SENSOR CONNECTION:** Wrap the sensor cord around the strain relief, and then connect to the "SENSOR PAD" receptacle on monitor.



- 5** **INCREASE VOLUME:** Increase volume by pressing and holding "VOL" button (solid green light indicates maximum volume).




- DECREASE VOLUME:** To lower volume press "VOL" button, release, press again and hold (solid yellow light indicates minimum volume)

Use



- 1** **TEST THE ALARM:** Put weight on sensor pad until you hear a double beep and see the green in-use light illuminate. Release weight to activate alarm.



- 2** **SILENCE THE ALARM:** return weight to the sensor pad or press "RESET".



- 3** **SECURE THE MONITOR** to the back of the patient's chair or bed, and (if applicable) connect the AC power cord and plug it into the room power receptacle.



- 4** **POSITION THE SENSOR PAD** on the chair or bed as instructed/illustrated on the pad.



- 5** **CONNECTING TO NURSE CALL:** Plug the nurse call connection cable into the "NURSE CALL" socket on the monitor, connect the opposite end of the cable into the Nurse Call Station (an adaptor may be required) and test.*



- 6** **CONNECTING A CALL CORD:** Plug the call cord into the "CALL CORD" jack in the monitor and test.*

* **NOTE:** NEVER connect a Nurse Call Connection Cable from the "NURSE CALL" socket on the side of the monitor to the "CALL CORD" socket on the other side of the monitor. Doing so will prematurely drain batteries.

©2015 Curbell Medical Products, Inc. MAP1588A

CSM-BC201 TROUBLESHOOTING CHECKLIST



This checklist has been developed to help identify common issues that can be quickly addressed to prevent unnecessary product returns. Please also refer to the product's Directions for Use and instructional video on Curbell's website (www.curbellmedical.com). If, after completing the checklist, the product continues to malfunction, please call your Curbell Representative at 888-222-2188.



No Power/Green "In Use" light does not blink

- Batteries are installed according to the diagram in the battery compartment.
- Batteries are fresh and/or known to be good.

If an AC adapter being used:

- The connection between the monitor and the AC outlet is secure and not damaged.
- The outlet in use is functioning properly and not controlled by a light switch.

Yellow "LOW BATTERY/REPLACE PAD" LED flashes

- The breakaway cord connector has been pushed in fully until it clicks.
- Batteries are fresh and/or known to be good.
- The pad/mat has been replaced.

No alert is given to the facility's nurse call system

- The connection between the monitor and the nurse call patient station is securely connected and not damaged (a 1/4" mono nurse call jack and/or an adaptor for compatibility is required).

False alarm or no alarm sound occurs

- The pad is placed properly on the chair or bed.
- No part of the pad/mat (including the neck of the pad) is folded and/or creased.
- The pad/mat has not expired.
- There is no heavy blanket, gel, or foam cushion on the pad that may interfere with proper operation or cause delay in the alarm.
- The volume level has been turned up.
- The Pause switch located in the battery compartment is not turned on. If it's in the ON position, when the reset button is pressed, the monitor is programmed to not alarm for 30 seconds.
- The delay setting (located in the battery compartment) is set at zero seconds.
- The cord between the monitor and sensor is securely connected and not damaged.
- The breakaway cord connector has been pushed in fully until it clicks.