


Gen4™ Pillow Speakers for Responder® 5



Enhanced (DirectAccess®)

Standard (Digital)



When we started designing the new Gen4 digital pillow speaker line, we wanted it to not only be a new platform for the future, but a new way of thinking, as well. Using over 45 years of experience designing and manufacturing the most popular pillow speakers in the world as our guide, we started from the ground up, combining the best features of our previous pillow speakers, suggestions from customers, and new materials and techniques to make Gen4 a step ahead of anything else in the nurse call industry today.

PATIENT SATISFACTION

- ▶ Buttons are easy to find and press.
- ▶ Streamlined and lightweight, like home remote controls.
- ▶ Crystal clear sound.

INFECTION CONTROL

- ▶ Proprietary design reduces the potential risk of infectious fluids or solids settling into the unit where they could be passed from patient to patient.
- ▶ Smooth shape allows for quick and easy cleaning, helping to ensure that dirt and grime cannot collect on the unit.

OXYGEN SAFETY

- ▶ UL Listed pillow speakers are tested with the nurse call system and verified for safe use in an oxygen-enriched atmosphere.
- ▶ Visit www.oxygensafe.com for more information

REGULATORY COMPLIANCE

According to 21 CFR 890.3710, powered communication systems and their accessories are considered Class II medical devices. As such, Curbell manufactures our nurse call communication accessories in compliance with this regulation. The FDA regulates all products and manufacturers within this regulation. Curbell's quality system complies with current Good Manufacturing Practices (cGMP) per 21 CFR 820 and is ISO 13485:2003 Medical Devices – Quality Management Systems certified.

FULLY SEALED UNIT

The sealed design reduces the potential risk of infectious fluids or solids settling into the unit.

AUXILIARY FUNCTIONS

Three auxiliary nurse call functions are available on all models. Up to four* environmental auxiliaries are available on Enhanced models, and two environmental auxiliaries are available on Standard models.

LARGE BUTTONS

Gen4's large buttons are recognizable at a glance and easy to press. The simplified layout makes contacting the Nurse and controlling TV functions simple.

USER-FRIENDLY BUTTON LAYOUT

Intuitive, standardized button positioning makes operation easy.

REPAIRABLE

You can utilize Curbell's repair service or choose to make basic repairs yourself, including quick cable replacement.

FLEX AND STRAIN RELIEF

Gen4's long flex relief was designed for longer life and greater flexibility. The smooth design facilitates easy cleaning and resists dirt and grime buildup. Internal strain relief prevents damage to internal components.

SUPERIOR SOUND

Clear, high quality sound for increased patient satisfaction.

OPTIONAL STEREO HEADPHONE JACK

Ideal for creating a quieter environment for patient rooms.



EASY-TO-CLEAN CASE

The streamlined design facilitates quick and easy cleaning.

VOLUME CONTROL OPTIONS

All models come standard with digital volume controls, and are available with analog volume for situations where two pillow speakers need to control one TV.



COMFORT-GRIP CASE

Gen4's lightweight case was designed to mimic the comfort and familiarity of a home remote.

MICROPHONE

The built-in microphone allows duplex audio.

	Part Number	Room Auxiliaries	Analog Volume	Digital Volume
Standard (DPS)*	350220	0	4D410U-001000A	4D410U-101001A
	350221	1 Light	4D410U-001040A	4D410U-101041A
	350222	2 Lights	4D410U-001030A	4D410U-101031A
Enhanced (DirectAccess)*	350200	0	4C410U-001000A	4C410U-101001A
	350201	1 Light	4C410U-001040A	4C410U-101041A
	350202	2 Lights	4C410U-001030A	4C410U-101031A

* Verify that your system is capable of supporting four auxiliaries