

TotalCare® and TotalCare SpO₂RT® Replacement Mattress Covers

Avoid the extra cost of replacing the complete mattress if the air bladders and foam of the original mattress are still acceptable — simply replace the worn-out cover with a new mattress cover from Curbell. These covers feature fluid-proof, waterproof welded seams, anti-microbial protection, and a waterfall flap over the zipper to keep just about anything from ruining your mattresses.

- ▶ Top cover:
 - Fade-resistant, polyurethane coated knit
 - Shear-resistant, breathable, bi-directional stretch urethane
 - Anti-microbial and hypo-allergenic fabric
 - BravoSilver Technology (SpO₂RT only) to help limit cross contamination (controlled release of silver ions resulting in a >99.9% reduction of bacteria)
- ▶ Tri-laminate vinyl bottom (TotalCare only) provides extra durable surface to withstand wear against bed frame
- ▶ Two-piece construction with full perimeter zipper for easier and faster replacement (can reduce replacement time by 50% or more), and easy maintenance (cleaning, changing bladders and hoses)
- ▶ Reinforced seams at high stress locations
- ▶ Bottom holes (with removable plugs) for air hoses and electronics (TotalCare: 11 holes, TotalCare SpO₂RT: 14 holes (classic pulmonary surface))
- ▶ Fully sealed cover prevents fluid ingress
- ▶ Easy to clean, stain-resistant, and durable
- ▶ Patented zipper-flap design to protect the mattress core
- ▶ Added fire retardants (meets CAL 117 fire rating)
- ▶ Five year warranty
- ▶ Also available: kits, foam topper, moisture/vapor barrier, fire barrier



Specially sized reinforced pneumatics exit ports (includes removable plugs)



Waterproof welded seams to help prevent cross contamination and fluid ingress



TotalCare

9257	36" x 84"	Urethane Cover, full four-side zipper	\$226.03
9259	36" x 84"	Kit (urethane cover, moisture/vapor barrier, foam topper)	\$398.00
9258	36" x 84" x 2"	Foam Topper, polyurethane	\$145.83
9255	36" x 84"	Moisture/Vapor Barrier	\$27.20
14648	36" x 84"	Fire Barrier	\$69.63

TotalCare SpO₂RT

9254	36" x 84"	Urethane Cover, full four-side zipper	\$576.72
9256	36" x 84"	Refresher Kit (Urethane Cover, Moisture/Vapor Barrier, Foam Topper, Fire Barrier)	\$648.80
9255	36" x 84"	Moisture/Vapor Barrier	\$27.20
14648	36" x 84"	Fire Barrier	\$69.63