

UNIVERSAL STRETCHER MATTRESSES AND COVERS

Curbell's stretcher mattresses are an economical replacement that can save you up to 50%. They feature high density foam, anti-microbial fabric, and are available in sizes that fit most stretcher frames. The tops are available in three fabrics:

POLYCARBONATE (blue)

- Non-conductive
- Resistant to most hospital disinfectants, including bleach and peroxide (disinfection data and a list of cleaners is available)
- Can be sterilized in an autoclave (134° C for 50 cycles)
- Our most durable option

POLYURETHANE (blue)


- Non-conductive
- Available with a viscoelastic ("memory foam") topper

LECTROLITE® (black)


- Conductive
- Helps control electrostatic discharge
- Prevents the accumulation of a static charge that would interfere with sensitive electronic equipment
- Thicker material

COMMON FEATURES


- Hypo-allergenic
- Anti-microbial
- Fade-resistant
- Shear-resistant breathable bi-directional stretch
- Replacement covers are also available separately
- Available with 0, 2, or 4 tapered corners
- 100% High density, anti-microbial, medical-grade polyurethane foam
- Non-slip PVC-coated poly knit bottom
- Sewn-in hook and loop positioners on the bottom
- Easy to clean, stain-resistant, and durable
- Fully sealed to prevent fluid ingress
- Zipper along one side (flap protects the foam core)
- Cover and foam have added fire retardants (meets CAL 117)
- Three year warranty




Mattresses are available with zero, two, or four tapered corners (four tapers shown)



Welded seams prevent fluid ingress



Non-slip bottom with sewn-in hook and loop positioner



Hook and loop positioner placement

SIZE	TAPERED CORNERS	POLYCARBONATE (NON-CONDUCTIVE)				POLYURETHANE (NON-CONDUCTIVE)				LECTROLITE® (CONDUCTIVE)					
		MATTRESS AND COVER		COVER ONLY		MATTRESS AND COVER		MATTRESS WITH FOAM TOPPER AND COVER		COVER ONLY		MATTRESS AND COVER		COVER ONLY	
3" x 26" x 76"	0					14576	\$245.02			16931	\$143.37				
	2	20736	\$266.00	20740	\$160.03	13597	\$245.02			16932	\$143.37	15323	\$254.71	16945	\$124.55
	4											14930	\$254.71	16946	\$124.55
3" x 29" x 76"	2					13598	\$254.36			16933	\$141.78				
4" x 24" x 76"	2	20737	\$285.92	20741	\$160.03	13599	\$262.58	14616	\$320.97	16934	\$143.37				
	0					14577	\$271.23			16935	\$141.78				
4" x 26" x 76"	2	19928	\$295.63	20742	\$160.07	13601	\$271.23	14617	\$360.36	16936	\$141.78	13600	\$263.19	16947	\$130.21
	4	20738	\$295.63	20743	\$160.07	13602	\$271.23			16937	\$141.78				
	0					13604	\$283.34			16938	\$141.78				
4" x 29" x 76"	2	19578	\$309.30	20744	\$160.03	13605	\$283.34	14619	\$381.53	16939	\$141.78	13603	\$271.66	16948	\$133.05
	4					13606	\$283.34			16940	\$141.78				
4" x 31" x 76"	2							14721	\$364.28	16941	\$141.78				
	4					18388	\$287.90								
5" x 26" x 76"	2	19579	\$326.70	19880	\$162.37	13607	\$323.63			16942	\$141.78				
5" x 29" x 76"	2	20739	\$336.93	20745	\$162.37	13608	\$322.59			16943	\$134.55				
	4					13609	\$322.59			16944	\$134.55				

Most common sizes and styles are listed here. Other sizes and styles are available.

Lectrolite is a registered trademark of Herculite, Inc.

©2016 Curbell Medical Products, Inc. MAP1351G