

# REPLACEMENT MATTRESS FOR STRYKER® GYNNIE® OB/GYN STRETCHER

Curbell's stretcher mattresses are an economical replacement that can save you up to 50%. They feature high density foam and anti-microbial fabric. The models listed below are specially designed to fit Stryker's Gynnie OB/GYN stretcher.

## FABRIC OPTIONS

### POLYURETHANE



- Non-conductive

### POLYCARBONATE

- Non-conductive
- Resistant to most hospital disinfectants, including bleach and peroxide (disinfection data and a list of cleaners is available)
- Can be sterilized in an autoclave (134° C for 50 cycles)
- Our most durable option

## COMMON FEATURES

- 4" thick comfort mattress
- 100% high density anti-microbial medical-grade polyurethane foam
- Cover and foam have added fire retardants (meets CAL 117)
- Hypo-allergenic
- Easy to clean, stain-resistant, and durable
- Fade-resistant
- Shear-resistant breathable bi-directional stretch
- Fully sealed to prevent fluid ingress
- Hook and loop placement and snap hardware are included
- Five year warranty



Fully sealed with waterproof welding to prevent fluid ingress



Attaches to stretcher deck with the included hook and loop fasteners and snaps

PART NUMBER	FABRIC	DESCRIPTION	PRICE
14610	Polyurethane	Stryker Gynnie Main Pad and Thigh Pad replacement set	\$444.97
14333	Polyurethane	Stryker Gynnie Main Pad section replacement (26"W x 59"L x 4"H/thick)	\$316.65
14334	Polyurethane	Stryker Gynnie Thigh Pad section replacement (26"W x 16.5"L x 4"H/thick)	\$128.32
19811	Polycarbonate	Stryker Gynnie Main Pad and Thigh Pad replacement set	\$534.45
21762	Polycarbonate	Stryker Gynnie Main Pad section replacement (26"W x 59"L x 4"H/thick)	\$376.79
21763	Polycarbonate	Stryker Gynnie Thigh Pad section replacement (26"W x 16.5"L x 4"H/thick)	\$157.66