

SAFE MOBILITY LIGHT

Insufficient lighting in a senior's room can lead to falls at the bedside at night.¹


Curbell Medical's Safe Mobility Light is an innovative solution that automatically activates lights when a senior leaves the bed.




The Safe Mobility Light offers tremendous safety benefits for individuals in assisted living or long-term care facilities. It's a cost-effective upgrade that can easily integrate into existing beds without decreasing comfort.



THE LIGHT BARS




The two light bars offer just enough illumination for safer bed exits. They attach magnetically or with included zip ties and can be installed anywhere on the bed frame.



THE SENSOR PAD



The sensor pad fits beneath the mattress, rather than under the sheet — increasing comfort, and eliminating the need for frequent cleaning.



EXPANDING THE SYSTEM



To expand coverage, the system can be enhanced with the addition of the BC600 Cordless Fall Monitor. Used in conjunction with additional sensors, this will extend wireless monitoring to the bed and other fall risk areas, such as a chair or the bathroom.

1. López-Soto P.J., Smolensky M.H., Sackett-Lundeen L.L., De Giorgi A., Rodríguez-Borrego M.A., Manfredini R., Pelati C., Fabbian F. Temporal patterns of in-hospital falls of elderly patients. Nurs. Res. 2016;65:435–445. doi: 10.1097/NNR.0000000000000184.