

CORDED BED AND CHAIR MONITOR

DESIGNED WITH CLINICIAN INPUT

The BC300 monitor is made for easy common day-to-day caregiver use. It facilitates timely mobility assistance to help keep your patients safe from falls by alerting staff that a patient has moved off of a sensor pad, which can indicate their intention to exit the bed or chair.

EASE OF USE

- Dedicated controls for Reset, Volume, Delay, and Tone
- “REPLACE BATTERY”, “SIGNAL LOSS”, and “REPLACE SENSOR” indicator LEDs show status at a glance

FLEXIBILITY

- Nurse call jack allows direct notification to nurse call
- Volume has three different audible levels with a silent alert option (when connected to nurse call) to reduce alarm fatigue and not cause the patient to restrain their movement (in compliance with CMS Guidance §483 11-28-17)*
- Delay control helps prevent false alarms when the patient momentarily moves off the sensor pad
- Tone switch allows to designate one of three settings according to facility preference or to dedicate to a specific patient
- AC adaptor option extends battery life, reducing maintenance and replacement costs
- Plastic strap and plastic wall mount bracket are included with the monitor (additional bed and chair bracket options are also available)

DURABILITY AND SAFETY

- Strain relief on the bottom of the monitor helps prevent damage to the cord when using with a corded sensor pad
- Durable thick plastic casing and impact-resistant protective rubber boot reduces replacement/repair costs
- Conforms to UL 1069 to ensure nurse call compatibility, safety, and durability
- Two year warranty



Volume, Tone, and Delay controls are kept out of sight in the battery compartment



Nurse call jack allows direct notification to nurse call



Cleanable plastic strap

* CMS Guidance §483 relates specifically to audible alarms in long-term care facilities