

















































ADAPTORS

The needs of some residents may require that multiple devices are connected to your nurse call system. Our line of adaptors will give your nurse call system the flexibility to handle this increased demand. They're designed to allow nurse call devices and Curbell's fall management products to connect to your nurse call system simultaneously.

- ▶ Multiple receptacle options
- ▶ Cable Strain relief
- ▶ Molded connectors
- ▶ Cord out notification
- ▶ Various styles, receptacles, and mounting flange options are available



Nurse Call System and Plug Style		Receptacle 1	Receptacle 2	Style	Part Number
	Executone Care/Com	 7-pin	 1/4" jack	Box	3248
	Executone Futura	 7-pin	 1/4" jack	Box	3254
	GE/Dukane System 2000, 2010, 2030, 2070	 13 Pin	 1/4" jack	Box	3271
	GE/Dukane ProCare 4000, 6000, System 2070, Telligence	 18-pin	 1/4" jack	Box	3363
	Hill-Rom COMposer, COMLinX, NaviCare	 20-pin	 1/4" jack	Box	3367
	Intego and TekTone	 8-pin Modular	 1/4" jack	Box	3322
	Jeron Provider 680	 8-pin DIN	 1/4" jack	Box	3283
		 8-pin DIN	 1/4" jack	Y Cable (with dummy plugs)	18194
	Rauland Responder III	 12-pin	 1/4" jack	Box	3266
	Rauland Responder III+, IV, 4000	 8-pin DIN	 1/4" jack	Box	3294
		 8-pin DIN	 1/4" jack	Y Cable (with dummy plugs)	18193
	SimplexGrinnell EZCare	 15-pin	 1/4" Jack	Box	3349
	Sylvania/HCE	 20-pin	 1/4" Jack	Box	3260
	West-Com	 10-pin Modular	 1/4" jack	Box	3341
	Universal 1/4"	 1/4" jack	 1/4" jack	Y Cable (with dummy plugs)	8507
		 1/4" jack	 1/4" jack	Y Cable	8506
	Universal 37-pin	 37-pin D-sub	 1/4" jack	Box	16659

NOTE: adaptor wiring requirements for nurse call and bed exit may vary. Please call us to discuss your application.