

PILLOW SPEAKER BREAKAWAY CABLE FOR RAULAND RESPONDER® 5



When a pillow speaker cable gets pulled from the patient station, there's potential for damage, which can cause several issues:

- ▶ Intermittent nurse call functionality
- ▶ False alarms
- ▶ Nurse call patient station repair costs
- ▶ An unavailable room during the repair

Our breakaway cable secures to the wall and safely breaks apart at a point away from the patient station, reducing cable replacement and patient station repairs.

- ▶ 9-Pin DIN Male to 9-Pin DIN Female
- ▶ Authorized by Rauland for their Responder 5 system
- ▶ Guide markings to match system and pillow speaker cable
- ▶ Cable tether included
- ▶ 2 foot length
- ▶ Extends overall cable length
- ▶ Part number JUP-9000-2

