

PEDIATRIC GEN4® PILLOW SPEAKERS

When your patients are children, special care needs to be taken to make sure that the environment they're in isn't intimidating. Curbell's Pediatric Gen4 Pillow Speakers help to bring a little fun to the patient's room, and make it a little more like home. With a bright, colorful design, and large buttons that make things like watching the TV and calling the nurse easy, the Pediatric Gen4 Pillow Speaker is a little thing that can make a big impact on your little patients and their families.

PATIENT SATISFACTION

Gen4 provides a quality TV experience that will help keep your pediatric patients satisfied.

- ▶ The buttons on the fun, kid-friendly overlay are easy to find and press
- ▶ Streamlined and lightweight case, like home remote controls
- ▶ Clean, clear sound

INFECTION CONTROL

Hospital-acquired infections can be very costly, so Gen4 is fully sealed¹ to help prevent infectious contaminants from getting into the case where they could be passed from patient to patient.

- ▶ Ingress Protection Rating of IP54 as defined by IEC 60529.*
- ▶ Smooth shape allows for quick and easy cleaning, and helps to prevent contaminants from collecting on the unit.

SOUND OPTIONS

- ▶ Optional headphone jack to help minimize noise
- ▶ Optional analog volume for situations where two devices (pillow speakers, beds, or Curbell Remote Entertainment Stations) control one TV.

EASE OF USE AND DURABILITY

- ▶ Simple reprogramming for different hospital-grade televisions
- ▶ High-impact case material and flexible cable.
- ▶ User-replaceable cable, speaker, headphone jack, and overlay

PATIENT ROOM INTEGRATION

The pillow speaker gives patients access to their care team, entertainment, and environmental controls, so we work closely with the manufacturers of these systems to make sure everything works together seamlessly.

- ▶ Up to three nurse call auxiliary buttons
- ▶ Up to four environmental auxiliary buttons

REGULATORY COMPLIANCE

In the US, pillow speakers are classified as medical devices, and manufacturers must be registered with the FDA. Curbell is fully compliant with all government and industry regulations. Gen4 is designed, manufactured, and tested to strict standards (like UL 1069) to give your patients a reliable connection, and to give you one less thing to worry about.



* For units without a headphone or keyboard jack

FULLY SEALED UNIT

The sealed design reduces the potential risk of infectious fluids or solids settling into the unit — even in the speaker area.

SUPERIOR SOUND

Clear, high quality sound to increase patient satisfaction.

AUXILIARY FUNCTIONS

The Gen4 platform allows for up to six auxiliary buttons to control lighting, temperature, drapes, nurse call requests, etc.

LARGE BUTTONS

Large buttons are recognizable at a glance and easy to press. The fun, kid-friendly layout makes contacting the Nurse and controlling TV functions simple.

DURABILITY

The high-impact case material is designed to withstand over 100 drops onto a concrete floor without impacting its functionality.

COMFORT-GRIP CASE

The lightweight ergonomic design mimics the comfort and familiarity of a home TV remote.

EASY-TO-CLEAN CASE

The streamlined design facilitates quick and easy cleaning with most hospital cleaning agents.

FLEX AND STRAIN RELIEF

Gen4's long flex relief was designed for longer life and greater flexibility. Internal strain relief prevents damage to internal components.

REPAIRABILITY

You can utilize Curbell's repair service or choose to make basic repairs yourself, including quick cable replacement.

ADDITIONAL FEATURES



OPTIONAL HEADPHONE JACK
Plug in 1/8" plug earphones to create a quieter environment.



ANALOG VOLUME CONTROL
Analog volume is available to ensure independent patient control in situations where two devices control one TV.



REPAIRABILITY
Gen4 is built to last; but if necessary, you can easily replace the cable, speaker, headphone jack, and overlay.

for more information, call your Curbell representative at **1-800-235-7500**