



Gen4
PILLOW SPEAKERS



YOU CAN COUNT ON CURBELL

FOR INNOVATION

From our first pillow speaker in 1962, Curbell has consistently led the industry with innovative features like nurse call integration, smart TV control, and extensive room controls. With easy-to-press buttons and clean design, Gen4+ Pillow Speakers are the latest in our long line of industry-leading patient control devices.



FOR EXPERTISE

Leading manufacturers of nurse call systems, healthcare TVs, and patient engagement systems depend on Curbell for their specialized private-labeled pillow speakers. When it comes to replacing pillow speakers in your hospital, count on Curbell to share the latest market insights and understand your specific needs to recommend the right unit for your patients.



FOR QUALITY AND SAFETY

In the United States, pillow speakers are classified as Class 2 exempt medical devices, and manufacturers must be registered with the FDA. Curbell is fully compliant with all government and industry regulations.

Gen4+ is designed, manufactured, and tested to strict standards, including UL 1069 for Hospital Signaling and Nurse Call Equipment, which includes intensive flammability, environmental, static discharge, and drop testing to ensure safe use for patients in hospital room settings.

In addition, Curbell ISO certified to 13485:2016. This standard specifies requirements for a quality management system where an organization needs to demonstrate its ability to provide medical devices and related services that consistently meet customer and applicable regulatory requirements.

EASY TO USE FOR PATIENTS OF ALL AGES

Knowing that patients aren't at their best either mentally or physically when they're in the hospital, Gen4+ Pillow Speakers are designed to be intuitive and simple to use.

VISIBILITY

The high-contrast case color and graphics make it easy for patients to find and use.

SPECIFIC NURSE CALL BUTTONS

Enable patients to relay their specific needs to the nursing staff.

DURABILITY

The high-impact case material is designed to withstand over 100 drops onto a concrete floor without impacting its functionality.

SEALED CASE

The sealed design reduces the potential risk of infectious fluids or solids settling into the unit. No exposed screws means one less area to attract dirt.

EASIER-TO-PRESS BUTTONS

Large buttons are recognizable at a glance and easy to press.

**2 YEAR
WARRANTY**

CRISP, CLEAR SOUND

Clear, high quality sound to increase patient satisfaction.

SOUND OPTIONS

Add a headphone jack to create a quieter environment and analog volume to ensure independent control in situations where two devices control one TV.



ROOM CONTROLS

The Gen4+ platform allows for up to four buttons to control lighting, temperature, shades, etc.

DASH BUTTON

Allows patients to access digital subchannels

EASY TO CLEAN

The streamlined design facilitates quick and easy cleaning to aid in infection control.


FLEX AND STRAIN RELIEF

The long flex relief is designed for longer life and greater flexibility. Internal strain relief prevents damage to internal components.

GEN4+ GIVES YOU OPTIONS

When your hospital has specific requirements for pillow speakers, you can count on Curbell to customize them to get you exactly what you need.

With over 50 years of experience and literally thousands of different models to choose from, we'll work with you to meet all your unique needs to help you deliver the best patient experience possible.



Gen4+ can even be customized to promote your hospital's brand.

PUT YOUR PATIENTS IN CONTROL

Knowing that patients aren't always at full strength when in the hospital, we've concentrated on making Gen4+ easy to use to give them a better experience and more control of their environment with less effort.

▶ NURSE CALL

Whether they're in pain, want something to drink, or need to use the bathroom, Gen4+ makes calling the nurse easy.

▶ THE ROOM ENVIRONMENT

Gen4+ can be customized to work with the systems in your rooms to control lighting, temperature, window shades, and more.

▶ TELEVISION

With a clean layout and clear sound, Gen4+ makes it easy for patients to pass the time with their favorite TV show.




BEDSIDE HOLDER
The bedside holder fits over most side rails and puts the pillow speaker within easy reach of the patient.



MINIMIZE REPAIR COSTS WITH PILLOW SPEAKER BREAKAWAY CABLES

A small investment designed to pay for itself in less than the cost of a single pillow speaker or patient station repair.

Pillow speakers need to hold up to the rigors of constant use by patients. To help reduce plug and patient station damage from excessive pulling, simply add a low cost Pillow Speaker Breakaway Cable which protects both by breaking apart cleanly when pulled in any direction.



REPAIR OPTIONS

CURBELL REPAIR SERVICE

Full testing and repair by Curbell technicians backed by a 90 day warranty.

DO-IT-YOURSELF

If you prefer to repair yourself, simply contact your sales representative to confirm parts required.


GEN4+ PILLOW SPEAKERS WITH BACKLIGHTING

Our Gen4+ pillow speakers for Rauland® Responder® 5 systems are now available with backlit buttons. Just plug the pillow speaker into the patient station and all the buttons will be clearly and evenly lit from behind to make them visible in any light — without negatively affecting nurse call system functionality. With a backlit Gen4+ Pillow Speaker, your patients will no longer be left in the dark — they'll have just enough light to see the buttons without keeping them up at night.





LEADING THE (R)EVOLUTION OF PATIENT CONTROL DEVICES



In hospital rooms around the world, patients are dealing with demands on their health, and staff members are dealing with demands on their time. Curbell Medical is at the center of it all, with product lines that have set the standard for hospital room integration. Put Curbell's experience to work for you to ensure the best possible patient experience.

RÊGO[®]
PATIENT INTERACTION SYSTEM



As patients desire greater access to clinical information, social connectivity, and additional entertainment options, many hospitals are providing commercial tablets to patients to enhance the hospital stay. These consumer-grade tablets introduce challenges with connectivity, charging, device tracking, and cleaning, as well as privacy and security.

As the leader in patient communication devices for over 50 years, Curbell Medical proudly introduces the new Rêgo™ Patient Interaction System – a scalable platform that addresses all the challenges associated with consumer-grade tablets, integrating all patient controls into a single, medical-grade device.

**For more information, contact your
Curbell representative, or visit
www.hellorego.com**

CURBELL
MEDICAL

www.curbellmedical.com • 1-800-235-7500