

GEN4+® PILLOW SPEAKERS

PUT YOUR PATIENTS IN CONTROL

Gen4+ Pillow Speakers give you the ability to give your patients a better experience in their room, and give them more control over a situation in which they might feel powerless. Knowing that patients aren't always at full strength when in the hospital, we've concentrated on making Gen4+ easy to use to give them more control of their environment with less effort.

▶ NURSE CALL


Whether they're in pain, want something to drink, or need to use the bathroom, Gen4+ makes calling the nurse easy.

▶ THE ROOM ENVIRONMENT

Gen4+ can be customized to work with the systems in your rooms to control lighting, temperature, window shades, and more.

▶ TELEVISION

With a clean layout and clear sound, Gen4+ makes it easy for patients to pass the time with their favorite TV show.



YOU CAN COUNT ON CURBELL

FOR INNOVATION

From our first pillow speaker in 1962, Curbell has consistently led the industry with innovative features like nurse call integration, direct access to TV channels, and extensive room controls. With easy-to-press buttons and clean design, our new Gen4+ Pillow Speakers are the latest in our long line of industry-leading patient control devices.

FOR EXPERTISE

Leading manufacturers of nurse call systems, healthcare TVs, and patient engagement systems depend on Curbell for their specialized private-labeled pillow speakers. When it comes to replacing pillow speakers in your hospital, count on Curbell to share the latest market insights and understand your specific needs to recommend the right unit for your patients.

FOR QUALITY

In the United States, pillow speakers are classified as Class 2 exempt medical devices, and manufacturers must be registered with the FDA. Curbell is fully compliant with all government and industry regulations. Gen4+ is designed, manufactured, and tested to strict standards, including UL 1069 for Hospital Signaling and Nurse Call Equipment, which includes intensive flammability, environmental, static discharge, and drop testing to ensure safe use for patients in hospital room settings.

* For units without a headphone or keyboard jack



©2019 Curbell Medical Products, Inc. MAPI754A

EASY TO USE FOR PATIENTS OF ALL AGES

Knowing that patients aren't at their best either mentally or physically when they're in the hospital, Gen4+ Pillow Speakers are designed to be intuitive and simple to use.

VISIBILITY

The high-contrast case color and graphics make it easy for patients to find and use.

SPECIFIC NURSE CALL BUTTONS

Enable patients to relay their specific needs to the nursing staff.

DURABILITY


The high-impact case material is designed to withstand over 100 drops onto a concrete floor without impacting its functionality.

FULLY SEALED

The sealed design reduces the potential risk of infectious fluids or solids settling into the unit — even in the speaker area.

EASIER-TO-PRESS BUTTONS

Large buttons are recognizable at a glance and easy to press.



CRISP, CLEAR SOUND

Clear, high quality sound to increase patient satisfaction.

SOUND OPTIONS

Add a headphone jack to create a quieter environment and analog volume to ensure independent control in situations where two devices control one TV.



ROOM CONTROLS

The Gen4+ platform allows for up to four buttons to control lighting, temperature, shades, etc.

DASH BUTTON

Allows patients to access digital subchannels

EASY TO CLEAN

The streamlined design facilitates quick and easy cleaning to aid in infection control.


FLEX AND STRAIN RELIEF

The long flex relief is designed for longer life and greater flexibility. Internal strain relief prevents damage to internal components.

GEN4+ GIVES YOU OPTIONS

When your hospital has specific requirements for pillow speakers, you can count on Curbell to customize them to get you exactly what you need.

With over 50 years of experience and literally thousands of different models to choose from, we'll work with you to meet all your unique needs to help you deliver the best patient experience possible.



Gen4+ can even be customized to promote your hospital's brand.

For more information, call your Curbell representative at **1-800-235-7500**