

Specialty Call Cords



Our Specialty Call Cords offer ease of use and comfort to those with special needs. They have ultra-sensitive touch surfaces that are ideal for patients with limited manual dexterity, and they can be used with most nurse call systems. Curbell offers five types of single and dual function Specialty Call Cords:



Breathcall®

- ▶ Pneumatic device activated by the patient's breath
- ▶ Single and dual models available
- ▶ For patients with special needs (severe burns, paraplegics, quadriplegics, immobilized patients)
- ▶ Adjustable clamp fits side rail or headboard
- ▶ Arm can be adjusted to almost any position
- ▶ Additional straws and filters available



E-Z Call™

- ▶ Can be placed anywhere for activation including the cheek, foot, hand, etc.
- ▶ Large tented design makes it easy for patient to see
- ▶ Linen clip securely attaches to bedding
- ▶ Velcro strip included for side rail mounting



Mechanical Pad

- ▶ Perfect for patients who may have difficulty using standard call cords or pillow speakers
- ▶ Large sensitive surface area
- ▶ Activates from any point on the pad's surface
- ▶ Linen clip securely attaches to bedding



Padcall®

- ▶ Large call surface (3" diameter)
- ▶ Sensitive, air-activated switch
- ▶ Gripper feet hold pendant in place
- ▶ Linen clip securely attaches to bedding
- ▶ Custom lengths available



Presscall®

- ▶ Air activated
- ▶ Fits comfortably in patient's palm
- ▶ Available in dual and single styles
- ▶ Linen clip securely attaches to bedding
- ▶ Custom lengths available

©2012 Curbell Medical Products, Inc. MAP292B