

Handheld Pendants

A Specialty Nurse Call Device

Curbell's handheld nurse call pendants are perfect for areas where a speaker isn't necessary, such as areas with swing-arm televisions, or emergency rooms and intensive care units that don't have televisions.

Ease of Operation

- ▶ Large, colorful buttons are easy for patients to understand and operate
- ▶ Raised texture on overlay and distinctive colors provide easy-to-locate buttons
- ▶ Buttons provide tactile feedback with an audible click
- ▶ Features a linen clip

Flexibility

- ▶ Designed to meet hospital's wiring requirements
- ▶ Custom overlays and cable lengths available
- ▶ Repairable
- ▶ Switch functions can be tailored to your hospital's needs
- ▶ Can support up to four auxiliary functions

Durability

- ▶ Long-lasting micro switches are rated for 30VDC, 1A
- ▶ Strain relief provided on all connections

Safety

- ▶ Safe for use in oxygen environments
- ▶ Electrostatic Discharge (ESD) resistant design
- ▶ Self-extinguishing case material



©2012 Curbell Medical Products, Inc. MAP589B