

Classic Call Cords



Our Classic Call Cords are designed to meet your facility's wiring requirements, and are safe for use in oxygen-enriched environments.

- ▶ Large contoured button is easy for patients to operate
- ▶ Button provides tactile feedback with an audible click
- ▶ Custom cable lengths available
- ▶ Switch rated for minimum of 30,000 cycles per UL 1069 section 24
- ▶ Electrostatic Discharge (ESD) resistant design
- ▶ Self-extinguishing pendant material



Call cords are available in numerous styles, including coiled, and dual.

©2012 Curbell Medical Products, Inc. MAP926B