

CABLES AND LEAD WIRES FOR THE GE® MAC® SYSTEM

When it comes to patient monitoring, hospitals don't want surprises. They want the reliability of Curbell Medical. A premier provider of patient monitoring solutions, our comprehensive line of high-quality cables and lead wires reflects Curbell's undisputed excellence in healthcare connectivity.

Curbell's ECG Cables, Lead Wires and Connector Clips are designed to meet or exceed the applicable performance requirements of ANSI/AAMI EC53A:/1998(R) 2001. Proper fit and function are assured, as each connector has been tested with OEM equipment to ensure compatibility and optimal performance.

- ▶ Fully shielded, low-noise cables minimize artifacts
- ▶ Quality ECG trace
- ▶ Standard AHA color-coded for quick and easy lead wire attachment
- ▶ ECG cables and lead wires are defibrillation-proof
- ▶ Designed in accordance with ANSI/AAMI EC53A: 2013 and ISO 13485:2016



| DESCRIPTION | PATIENT END | COLOR CODE | LENGTH | OEM PART NUMBER | CURBELL PART NUMBER |
|---------------------------------------------------------------------------------------------------------------------------|-------------|------------|---------|-----------------|---------------------|
|  Trunk Cable for GE CAM 14 | NA | AHA | 4' | 2016560-001 | 28080 |
| | | | 15' | 2016560-002 | 28081 |
|  GE MAC ECG Cable Trunk Cable, 10 Lead | NA | AHA | 8' | 2104723-001 | 28078 |
|  GE MAC ECG One-Piece Cable | Banana Clip | AHA | 8' | 2104727-001 | 28079 |
|  GE MAC ECG Lead Wire Set, Add-On, HEI | Banana Clip | AHA | 29" | 2104751-001 | 27467 |
|  GE MAC ECG Lead Wire Set, Base 10 | Pinch | AHA | 26"/46" | 2104749-001 | 27486 |
|  GE MAC ECG Lead Wire Set, Base 10 | Banana Clip | AHA | 26"/46" | 2104724-001 | 27419 |